

HEALTH HUB

Connects you with your care team

Plug it In

- » Plug in the Health Hub Console.

Step
01



Power On

- » Flip the Power Switch to the "ON" position (Located on the bottom of Console).
- » Your Health Hub will automatically connect to our cellular network.

Step
02



Press Call button

- » Press the lower left hand blue button to test the connection and make sure that it's properly working.

Step
03



WEIGHT SCALE

Measures body weight

Place the scale

- » Place scale on hard, flat surface. Do not use on carpet or uneven surfaces. Scale will turn on automatically.

Step
04



The results of your reading

- » Your Health Hub system will announce the results of your reading.

Step
05



Step off device

- » Step off device which will turn off automatically.

Step
06



**YOU'RE READY
AND PROTECTED!**



www.nsightcare.com